


The Audicus Aura CIC Hearing Aid Specifications

Fitting Range




Dimensions


length: 17.6 mm (0.7 in)

height: 11.1 mm (0.4 in)


width: 7.4 mm (0.3 in)



OSPL90 CURVE



FREQUENCY RESPONSE CURVE



ANSI, S3.22-2003, 2cc Coupler

OSPL90 Maximum	107 dB SPL
HF Average OSPL90	103 dB SPL
Peak Gain	32 dB
HFA Full-On-Gain	27 dB
Reference Test Gain	26 dB
Frequency Range	200 - 8000 Hz
THD at 500 Hz	<3%
800 Hz	<3%
1600 Hz	<3%
Equivalent Input Noise	<27 dB SPL
Attack Time @ 2kHz*	6-9 ms / 35-105 ms
Battery Current Drain	0.75 mA
Battery Life [10A]	105 hrs

