


# AUDICUS

## The Partner Mic for The Audicus Spirit Hearing Aid

### Description

- Directly stream audio from the Partner Mic to your Spirit hearing aids
- Allow crystal clear conversations in a noisy environment, or to hear a speaker from a distance
- Operates up to 25 meters / 80 feet
- Non-directional: attach anywhere close to the speaker with the convenient clip
- Connects to one pair of hearing aids at a time with an easy out-of-the-box pairing process
- Stream up to 6 hours with an internal lithium polymer battery
- Conveniently recharge with the included USB-C cable, ~2 hours
- RF Power: 20 mW



Height: 18mm (.7 in)  
Depth: 55mm (2.14 in)  
Width: 30mm (1.2 in)  
Weight: 18g (0.64oz)