


The Audicus Mini CIC Hearing Aid Specifications

Fitting Range




Dimensions

Length: 17.4 mm (0.7 in)

Height: 11.9 mm (0.5 in)


Width: 5.6 mm (0.2 in)




Max. Output sound pressure level

($L_1 = 90$ dB)

Full on gain
($L_1 = 50$ dB)



Frequency response ($L_1 = 60$ dB)



2cc Coupler

OSPL90 Maximum 114 dB SPL

HF Average OSPL90 108 dB SPL

Peak Gain 50 dB

HFA Full-On-Gain 45 dB

Reference Test Gain 31 dB

Frequency Range 100 - 8600 Hz

THD at 500 Hz <2%

800 Hz <3%

1600 Hz <2%

Equivalent Input Noise <20 dB SPL

Attack Time @ 2kHz* <15 ms

Battery Current Drain 1.2 mA

Battery Life [10A] ~70 hrs

